



Healing Biologix

Innovating Health Solutions

INTRODUCING HEALING BIOLOGIX: Your Physician-Supplied Wound Care System

Welcome to Healing Biologix. We are dedicated to providing innovative and convenient surgical dressings to patients like you. Your care provider has selected Healing Biologix products to maximize your healing potential, protect your incisions or wounds, and create a barrier against infection and other contaminants. Our expertise ensures that your surgical dressings are tailored to your post-operative wound care needs and are shipped directly to your home, arriving within 24 – 48 hours.

What to Expect: Your primary dressing will consist of Collagen Sheets or Collagen powder to cover the application area daily until the supplies are exhausted or your wound is healed.



PRIMARY DRESSING OPTIONS:

- **Collagen Sheets:** Our collagen dressings are porous membranes designed to be permeable and breathable, with fluid control similar to skin. They are flexible, pliable, elastic, and durable for easy application.
- **Collagen Powder:** The particles in our collagen powder have a high surface area, allowing them to effectively penetrate and deliver the maximum amount of collagen into the wound surface.

SECONDARY DRESSING OPTIONS:

- **Silicone Composite Dressing:** These absorbent dressings feature a non-traumatic border and provide an ideal cover for various wounds. Indicated as a secondary dressing to cover primary dressings, they assist in protecting the wound and preventing bacterial contamination or maceration. These dressings are also waterproof.
- **Island Composite Dressing:** These dressings feature a soft, absorbent pad that collects exudate and protects the wound from infection. The non-woven backing with an adhesive border is gentle and holds the pad in place. The island composite dressing also has a waterproof backing to keep wounds dry and remove external contaminants.

At Healing Biologix, we're committed to improving outcomes, compliance, and convenience for patients like you. If you have any questions about our products or their applications, please don't hesitate to contact us. Experience the difference with Healing Biologix – your partner in wound care.

To learn more, visit the Healing Biologix Website.



healingbiologix.com



Healing
Biologix

Innovating Health Solutions

COLLAGEN WOUND CARE DRESSING: INSTRUCTIONS

Postoperative Dressing Changes & Materials

You will receive a package from Healing Biologix containing supplies for changing your surgical incision dressings after surgery, typically arriving 2-3 days post-operation. These supplies are carefully selected to enhance your healing potential, protect and strengthen your incisions, and form a barrier against infection and contaminants.

INSTRUCTIONS FOR APPLICATION:

STEP 1: APPLY COLLAGEN DRESSING

For Longer Surgical Incisions (Collagen Sheet):

- Open the package and remove the collagen sheet.
- Cut the collagen sheet to fit and cover the incision, leaving approximately 2 inches of collagen border around all sides.
- Place the cut collagen sheet directly over the incision.

For Short or “Stab” Incisions (Collagen Powder):

- Open one packet of collagen powder for each incision.
- Apply all the collagen powder directly onto the incision.

STEP 2: APPLY SECONDARY DRESSING

After applying the collagen dressing, cover the area with either a “Silicone Composite Dressing” or a “Silicone Composite Island Dressing” to lock in the collagen and seal the wound.

Prior to Applying Collagen:

1. Completely remove the previous dressing(s) from the incision(s).
2. Cleanse the Incision(s) by using sterile saline or warm water with antibacterial soap to clean the incision(s). ***Avoid using peroxide unless necessary to remove dried blood.***
3. Dry the Incision(s) by gently patting the incision(s) dry with a clean towel, gauze, or use a hair dryer on a cool setting.

Additional Notes:

- The silicone dressings are water-resistant, allowing you to shower with the dressing from the previous day. Afterward, replace with a fresh set of collagen/silicone dressings.
- If you have any questions or concerns about the application of Healing Biologix Collagen Wound Care Dressings, please consult your healthcare provider or contact Healing Biologix directly.

Important:

- Follow these instructions carefully for optimal wound healing, and consult your healthcare provider if you experience any issues or concerns.

FOR MORE INFORMATION:

Call (844) 424-6564 | support@healingbiologix.com



healingbiologix.com